

# SYREBO® SY-UEA2

**Upper Limb Rehabilitation Robot**

**Advanced Gamified Rehabilitation  
for Arm and Shoulder Recovery**



# Revolutionizing Upper Extremity Rehabilitation


The Syrebo® Upper Limb Rehabilitation Robot SY-UEA2 combines flexible robotics, high-precision optical positioning technology, and engaging interactive games to provide effective, fun, and motivating training for arm and shoulder function recovery.

**Key Innovation:** By merging advanced robotics with neuroscience and gamification, the SY-UEA2 transforms repetitive therapy into an engaging experience that maximizes patient motivation, treatment adherence, and clinical outcomes.





# Core Features



## High-Precision Optical Positioning Technology

### Millimeter-Level Accuracy


- Advanced optical sensors track every movement
- Real-time monitoring of motion parameters
- Intelligent recognition of user status
- Precision reaching millimeter (mm) level
- Ensures accurate measurement and documentation



## Ergonomic Design

### Patient-Centered Engineering

- Adjustable arm length for wide range of patients
- Specialized accessory combination
- Compatible with wheelchairs
- One-touch height adjustment
- Suitable for multiple positions



## Omnidirectional Mobile Chassis

### Freedom of Movement

- Full-featured mobile platform
- Extensive training range
- Portable and compact design (only 8kg)
- Adaptable to multiple scenarios (sitting, standing, bedside)
- Easy movement between treatment areas



## Multidimensional Training System

### Comprehensive Rehabilitation

- 3 training goals: Motor, Cognition, ADL
- 4 motor training modes: Passive, Assistive, Active, Resistive
- Covers entire rehabilitation cycle
- Cognitive training integrated
- Activities of Daily Living (ADL) practice



# Motor Training 4 Progressive Modes

## Comprehensive Coverage of Rehabilitation Cycle

### ✓ **Passive Training (MMT Gr.0)** *For patients in early recovery stage*

- Robot fully guides arm movement
- No patient effort required
- Maintains range of motion
- Prevents contractures and stiffness

### ✓ **Assistance Training (MMT Gr.1-2)** *For patients with minimal movement*


- Robot provides partial assistance
- Patient contributes some effort
- Gradually increases patient participation
- Builds strength progressively

### ✓ **Active Training (MMT Gr.3)** *For patients with fair strength*

- Patient initiates and performs movement
- Robot follows the intended trajectory
- Provides guidance and support as needed
- Encourages maximum patient effort

### ✓ **Resistance Training (MMT Gr.4-5)** *For patients in late-stage recovery*

- Robot provides controlled resistance
- Strengthens muscles
- Challenges patient capabilities
- Prepares for functional independence



## Visual Feedback:

- Host trajectory (blue line): Robot's guided path
- Patient trajectory (orange line): Patient's actual movement
- Real-time comparison for immediate feedback

# Cognitive Training

## Enhance Neuroplasticity Through Engaging Activities

Improve user's cognitive function and enhance neuroplasticity through exercises targeting memory, attention, perception, and other functions.

### Attention Training "Catch Bees" Game

- Divided attention exercises
- Visual tracking practice
- Sustained attention building

### Sensory and Perceptual Training "Categorization" Game

- Colour recognition exercises
- Shape recognition challenges
- Visual discrimination practice

### Memory Training "Let Goods Back" Game

- Spatial memory development
- Focused attention training
- Working memory enhancement

### Unilateral Neglect Training "Something New" Game

- Working memory challenges
- Focused attention improvement
- Spatial awareness enhancement







# ADL Training

## Restore User's Ability for Daily Living

By practicing in real-life situations, users gradually improve their attention, strength, and ability to carry out daily activities.

### Functional Task Practice

#### Tidy the Room

- Attention training
- Muscle strength development
- Simulates organizing objects
- Real-world application

#### Shopping

- Attention enhancement
- Muscle strength building
- Grocery selection simulation
- Practical life skills

### Benefits of ADL Training

- Bridges therapy to real-world function
- Improves independence
- Increases patient engagement
- Provides meaningful context
- Prepares for home environment



# Diversified Customizable Training Programs




## Pre-Set Trajectories

### Joint-Specific Movements

- Protraction & retraction of scapula - Shoulder blade mobility
- Elbow joint flexion & extension - Bending and straightening
- Internal & external rotation of shoulder joint - Rotational movements

### Complex Movement Patterns

- Range of motion training - Full arc movements
- Unilateral neglect training - Spatial awareness
- Balance training - Coordination and stability



## Customize Trajectory

### Patient-Specific Programming

- Create unlimited custom movement patterns
- Adapt to individual patient needs
- Adjust complexity as patient progresses
- Combine multiple movement types
- Save personalized programs.



# Quantifying the Training Process

## Intelligent Data Analysis & Evaluation

### Comprehensive Progress Tracking

- Support multi-project training planning
- Combine evaluation with training data
- Quantitatively manage the rehabilitation process
- Generate detailed reports
- Evidence-based treatment decisions

## Multi-Dimensional Feedback



### Real-Time Monitoring

- Speed tracking during training
- Strength assessment
- Range of motion measurement
- Endurance evaluation
- Symmetry analysis

Therapist Benefits: According to real-time feedback, therapists can adjust the rehabilitation training plan in time to achieve closed management of the process.

### Patient Benefits:

- Visual progress tracking
- Clear goals and milestones
- Motivation through visible improvement
- Objective measurement of gains





# All-Round Safety Protection

## 4 Safety Systems




### Force Limitation Protection (E10)

- Prevents excessive force application
- Protects patient from strain
- "Please stop force overly to avoid injury" alert



### Overload Protection (E2)

- Monitors system capacity
- "Please restart the device or contact technical support" notification
- Prevents equipment damage



### Emergency Stop Function (E8)

- Activated during safety periods
- "Please stop training manually after safety period" warning
- Immediate halt capability



### Electronic Fence Touched

- Boundary protection system
- Prevents unintended movements
- Keeps training within safe zones

# Indications for Use

## Patient Populations

*Patients with upper limb dysfunction caused by:*



### Neurological Conditions

- Stroke (CVA) - Arm and shoulder recovery
- Brain injury - Traumatic brain injury rehabilitation
- Cerebral palsy - Pediatric and adult treatment
- Parkinson's disease - Movement and rigidity improvement
- Alzheimer's disease - Cognitive and motor maintenance



### Trauma and Injury

- Hand trauma - Arm and shoulder injuries
- Post-surgical rehabilitation - Orthopedic procedures
- Fractures - Upper extremity recovery



### Other Conditions

- Lymphedema - Swelling reduction, mobility improvement
- Muscle weakness - Deconditioning recovery
- Limited range of motion - Stiffness and contractures

# Suitable Clinical Environments

## Multiple Department Applications



### **Rehabilitation Department -**

Primary treatment setting



### **Neurosurgery Department -**

Post-operative rehabilitation



### **Neurology Department -**


Neurological condition management



**Geriatric Department -** Elderly patient care



**Pediatric Department -** Children's rehabilitation (pediatric version available)



**Orthopedic Department -** Musculoskeletal recovery



**Physical Therapy Department -** Outpatient rehabilitation



# Portable & Compact Design - Easy to Use

Meeting Various Scenarios: Adjustable for Multiple Positions



## ✓ Bedside Use

- Perfect for hospitalized patients
- Acute phase rehabilitation
- Intensive care settings
- Post-surgical early intervention

## ✓ Standing Position

- For patients with greater mobility
- Advanced rehabilitation stage
- Functional training emphasis
- Balance integration

## ✓ Sitting Position

- Standard outpatient clinic setup
- Most common configuration
- Comfortable for extended sessions
- Wheelchair compatible

## ✓ One-Touch Height Adjustment

- Compatible with barrier-free wheelchairs
- Smooth height transitions
- No manual lifting required
- Accessible for all patients



# Innovative Structure

## Omnidirectional Moving Chassis

### Training More Freely

- 360-degree movement capability
- Extensive workspace coverage
- Natural arm movement patterns
- No spatial limitations

## High-Precision Optical Positioning Technology

### Accuracy Reaching Millimeter (mm) Level

- Infrared sensors track position continuously
- Real-time coordinate mapping
- Precise trajectory following
- Microsecond response time

### Specialized Features

- Ergonomic design
- Specialized accessory combination
- Adjustable arm length
- Suitable for a wide range of patients
- 4 sets of safety protection systems
- Increasing safety assurance

# Evidence-Based Approach

## Scientific Foundation

Based on recent clinical studies, upper extremity neurorehabilitation effectively improves arm and hand function by leveraging:

### Key Principles:

- Intensive training - High-dose practice
- Task-specific activities - Meaningful movements
- High-repetition exercises - Neuroplasticity stimulation

**Emerging Technology Benefits:** Robotics, combined with virtual reality (VR) and engaging interfaces, show promise as effective complements or alternatives to conventional therapy, particularly for severely impaired patients.

### Clinical Evidence:

- Dose and intensity matter: Higher dose of task-oriented training is superior to lower-dose conventional approaches for improving upper extremity motor deficits
- Applicable across phases: Effective in both subacute and chronic phases after stroke
- Motivation increases adherence: Gamification significantly improves patient engagement and treatment compliance



# Why Choose Robotic Upper Extremity Rehabilitation?

## Advantages Over Traditional Therapy

### ✓ Consistent Dosage

- Robot provides exact same assistance/resistance every repetition
- Human therapists may fatigue or vary

### ✓ More Treatment Options

- Motor training (4 modes) + Cognitive training + ADL training
- Traditional therapy typically focuses on motor only

### ✓ Higher Repetition Count

- Patients can perform hundreds of repetitions per session
- Manual therapy limited by therapist endurance

### ✓ Increased Motivation

- Engaging games make therapy fun
- Traditional exercises can become monotonous

### ✓ Better Functional Tracking

- Detailed progress reports
- Visual proof of improvement

### ✓ Data-Driven Decisions

- Evidence-based treatment adjustments
- Clear documentation for insurance and outcomes

### ✓ Objective Measurement

- Precise tracking of every parameter
- Traditional therapy relies on subjective assessment



# Technical Specifications

## System Specifications

### More Treatment Options

Weight	8 kg (portable)
Chassis Type:	Omnidirectional mobile
Positioning Technology:	High-precision optical (millimeter leve)
Workspace:	Extensive range, all directions

### Training Capabilities

Training Goals:	3 (Motor, Cognition, ADL)
Motor Training Modes:	4 (Passive, Assistive, Active, Resistive)
Cognitive Games:	Multiple interactive options
ADL Scenarios:	Real-life task simulations
Pre-set Trajectories:	Multiple joint-specific patterns
Custom Trajectories:	Unlimited programmability

## Features

### More Treatment Options

Height Adjustment	Yes, compatible with wheelchairs
Use Positions:	Sitting, standing, bedside
Safety Systems:	4 comprehensive protection features
Real-time Monitoring:	Speed, strength, range, endurance
Report Generation:	Automatic with data analysis
Available Versions:	Adult and pediatric

### Regulatory

CE Certified:	European medical device approval
FDA Registered:	United States authorization
Warranty:	2 years standard





# Pediatric Version Available

## Syrebo® Rehabilitation for Children

### Specialized Design for Pediatric Patients

Syrebo has introduced a pediatric version of its rehabilitation robot, designed to assist children in hand and arm rehabilitation.

#### Ideal For:

- Pediatric cerebral palsy
- Traumatic brain injuries in children
- Congenital neurological conditions
- Post-surgical pediatric recovery
- Developmental delays

#### Child-Friendly Interface:

- Colorful, engaging games
- Age-appropriate challenges
- Shorter session durations
- Reward systems
- Parental involvement features

#### Special Features:

- Combines flexible robotics technology with neuroscience
- Helps with finger flexion and extension
- Reduces muscle tension
- Alleviates swelling and stiffness
- Promotes brain neuroplasticity
- Improves hand mobility
- Accelerates the rehabilitation process



# Clinical Settings

## Facilities Using the SY-UEA2

- Skilled Nursing Facilities (SNF) and Assisted Living (ALF) Use with elderly and frail persons requiring upper extremity rehabilitation
- Outpatient Neurological Training Centers Intensive therapy for stroke and brain injury patients
- Inpatient Rehabilitation Facilities Early intervention in acute and subacute phases
- Physical and Occupational Therapy Clinics Integration into comprehensive rehab programs
- Hospital Rehabilitation Departments Neurology, stroke, orthopedic departments
- Pediatric Rehabilitation Centers Specialized treatment for children with neurological conditions

# Immersive Interaction Experience

Engaging Gamified Therapy. The Syrebo SY-UEA2 uses colorful, interactive games to:

## Provide Immediate Visual Feedback

- See results in real-time
- Understand performance instantly

## Challenge Cognitively

- Train brain while training movement
- Dual-task activities for neuroplasticity

## Create Session Variety

- Multiple game options prevent boredom
- Progressive difficulty levels

## Measure Progress Objectively

- Game scores correlate with function
- Clear improvement metrics



## Example Game Scenarios:

- Underwater object collection
- Grocery shopping simulation
- Landscape navigation
- Object categorization
- Memory challenges
- And many more...



# Complete Rehabilitation Ecosystem



The SY-UEA2 is part of Syrebo's comprehensive rehabilitation solution portfolio:

## Hand Rehabilitation

- Robotic hand rehabilitation gloves
- Hand rehabilitation system with BCI (Brain-Computer Interface)
- Soft exoskeleton hand devices

## Upper Limb

- Upper Limb Rehabilitation Robot SY-UEA2
- Dynamic Arm Support systems
- Shoulder training devices

## Lower Limb

- Soft exoskeleton walkers
- Gait rehabilitation systems

## Complementary Therapy

- Low and medium frequency electrotherapy systems
- Nerve stimulation devices

## Global Presence & Trust

### Worldwide Impact

- 80+ Countries & Regions Syrebo technology available across six continents
- 4,000+ Medical Units Hospitals, clinics, and rehabilitation centers using Syrebo systems
- 60,000+ Families Patients and families benefiting from Syrebo rehabilitation

## Manufacturing Excellence

### *Shanghai Siyi Intelligent Technology Co., Ltd.*

- Founded with focus on intelligent rehabilitation robotics
- R&D team with 10+ years of experience in robotics, rehabilitation medicine, and AI
- Collaboration with world-renowned rehabilitation experts
- Commitment to effective, safe, and affordable solutions

### Quality Certifications:

- CE-MDR Certified
- FDA Registered
- ISO Standards Compliant

# Comprehensive Support

## by Neuro Rehab Recovery

*As an authorized distributor of Syrebo in the United States, Neuro Rehab Recovery provides:*

### Clinical Training

#### Expert Instruction

- Training by certified neuro/ortho trained therapists
- Comprehensive operational training
- Best practice protocols
- Clinical application strategies

### Ongoing Clinical Support

#### Continuous Access

- Neuro-trained physical therapist consultation
- Treatment optimization guidance
- Case discussion and problem-solving
- Protocol development assistance

### Technical Support

#### USA-Based Service

- Repair and technical support resources
- Prompt response times
- Replacement parts availability
- Remote troubleshooting



### Warranty

#### Quality Assurance

- 2-year standard warranty
- Reliable service and support
- Peace of mind protection

#### Flexible Business Models

- Purchase options
- Rental programs available
- Leasing arrangements
- Customized solutions for your facility



# Why the Syrebo® SY-UEA2?

## Complete Upper Extremity Solution

### High-Quality Manufacturing

- Precision engineering
- Durable components
- Reliable performance

### Fun and Easy to Use

- Intuitive interface
- Engaging for patients
- Simple for therapists

### Functional Progression Tracking

- Objective measurements
- Clear progress documentation
- Evidence for insurance

### Extensive Movement Options

- Multiple trajectory patterns
- Customizable programming
- All major arm movements

### Comprehensive Training Modalities

- Motor + Cognitive + ADL
- Full rehabilitation spectrum

### Adjustable for All Patients

- Height adjustment
- Wheelchair compatible
- Bedside, sitting, standing



### Pediatric and Adult Versions

- Solutions for all ages
- Specialized designs

### Strong Warranty and Support

- 2-year coverage
- USA-based service
- Clinical consultation

### Evidence-Based Approach

- Supported by research
- Proven clinical outcomes

# Elevate Your Upper Extremity Rehabilitation Program

The Syrebo® Upper Limb Rehabilitation Robot SY-UEA2 offers advanced technology, evidence-based training, and superior patient experience.

## Experience the Difference

**Schedule a Demonstration and ask about our flexible financing programs**

See the SY-UEA2 in action and discover how gamified robotic rehabilitation can transform your practice and your patients' outcomes.

***Empowering Lives Through Advanced Neuro Rehabilitation***



## Contact Information:

 **1 (866) 476-6897**

 [info@neurorehabrecovery.com](mailto:info@neurorehabrecovery.com)

 [www.neurorehabrecovery.com](http://www.neurorehabrecovery.com)

