

## Contact us

Email: info@neurorehabrecovery.com

Telf: 1 (866) 476 - 6897


Website: www.neurorehabrecovery.com













## MOTOmed gracile12 grows with you

The MOTOmed gracile12 grows with the child thanks to a wide range of adjustment options. For example, the pedal axle with foot shells is infinitely height adjustable. The MOTOmed gracile 12 is available as a leg trainer but also as a leg or arm/upper body trainer to meet the different patient requirements.

Non-binding fig. of MOTOmed gracile12 leg trainer

# MOTOmed<sub>®</sub> gracile12

### Uniquely child-friendly


The MOTOmed gracile12 was specially developed for children and small people with limited mobility. From a height of about 90 cm, children and adults can enjoy the efficient MOTOmed movement therapy with safety functions, motivational games and therapy programs. If required, the leg trainer can be extended with the arm/upper body trainer and other accessories. The quiet motor is conducive for both concentrated therapy sessions as well as for exercising at home.

#### Basic equipment


- · Pediatric safety foot shells and pediatric leg guides with calf shells
- · Smaller inner diameter of 12 cm
- · Sturdy full metal construction, high quality and stable
- · Large, user-friendly display (5.7")
- · Height adjustment of the pedal axle from 26 to 46 cm
- · 2-stage pedal radius adjustment (3.5 cm or 7 cm)
- · Easily movable with large castors (Ø 13 cm)
- Telescopic front stand
- · Colour: blue/white/chrome

#### Functions

- · Passive, Assistive and Active Training
- Electronic foot insertion aid
- Speed from 1 to 60 rpm (passive)
- · Resistance from 0 to 20 (active)
- · Motor power from 1 to 10 (passive)
- · Presetting of therapy time from 0 to 120 minutes
- MovementProtection & SpasmControl
- Quiet, gentle and harmonious cycle movement (SmoothDriveSystem)
- · Gentle starting and stopping of the movement
- Safety stop and SymmetryTraining
- Therapy and motivation programs
- · Detailed biofeedback during and after training
- · Games
- Language selection



The pedal axle of the leg trainer is steplessly height adjustable from 26 to 46 cm (measured from floor).



The child-friendly pedal distance of 12 cm avoids misalignment of the hip joint.

# MOTOmed<sub>®</sub> gracile12

### Uniquely child-friendly


The MOTOmed gracile12 was specially developed for children and small people with limited mobility. From a height of about 90 cm, children and adults can enjoy the efficient MOTOmed movement therapy with safety functions, motivational games and therapy programs. If required, the leg trainer can be extended with the arm/upper body trainer and other accessories. The quiet motor is conducive for both concentrated therapy sessions as well as for exercising at home.

#### Basic equipment


- Pediatric safety foot shells and pediatric leg guides with calf shells
- · Smaller inner diameter of 12 cm
- · Sturdy full metal construction, high quality and stable
- · Large, user-friendly display (5.7")
- · Height adjustment of the pedal axle from 26 to 46 cm
- · 2-stage pedal radius adjustment (3.5 cm or 7 cm)
- · Easily movable with large castors (Ø 13 cm)
- · Telescopic front stand
- · Colour: blue/white/chrome

#### Functions

- · Passive, Assistive and Active Training
- · Electronic foot insertion aid
- · Speed from 1 to 60 rpm (passive)
- · Resistance from 0 to 20 (active)
- · Motor power from 1 to 10 (passive)
- · Presetting of therapy time from 0 to 120 minutes
- MovementProtection & SpasmControl
- · Quiet, gentle and harmonious cycle movement (SmoothDriveSystem)
- · Gentle starting and stopping of the movement
- · Safety stop and SymmetryTraining
- Therapy and motivation programs
- · Detailed biofeedback during and after training
- · Games
- · Language selection



The pedal axle of the leg trainer is steplessly height adjustable from 26 to 46 cm (measured from floor).



The child-friendly pedal distance of 12 cm avoids misalignment of the hip joint.


# MOTOmed® gracile12

# Leg or arm/upper body trainer

### نے نے

Order number:

594.003 + 599.000



### Basic equipment (For all details, see equipment and functions page 23)

Pediatric safety foot shells (pair, incl. pedal cranks) Pediatric leg guides with calf shells (pair) Soft Grip handles (pair, incl. pedal cranks)

### Technical data

Dimensions (L×W×H):	60-88 × 55 × 99-114 cm
Weight:	33 kg
Maximum patient weight:	90 kg
Display size:	5.7" / 14.5 cm
Classification:	Ila according to MDD 93/42/EEC / MDR (EU) 2017/745

Mains voltage: 100-240 V~ / max. 120 VA, mains frequency: 47-63 Hz, class II equipment / Type BF

### Frequently used accessories (More accessories from page 40)

Designation	Order No.	Page
Plastic-coated pediatric safety foot shells (pair, incl. pedal cranks)	152.120	p. 44 (i)
Pediatric leg guides with plastic-coated calf shells (pair)	170.100	p. 45 🛈
Pedal radius quick adjustment for pediatric safety foot shells	588.020	p. 46 (i)
(leg trainer, pair)		
Self-operated foot holders (pair)	598.000	p. 45 (i)
Pediatric forearm shells with arm cuffs, size S (pair)	556.024	p. 52 (i)
Hand fixation by wrist cuff, size S (piece)	562.020	p. 48 🗓
Pediatric tetra handles (pair)	555.020	p. 50 (i)
Vertical handles with hand rests for children (pair)	557.020	p. 51 (i)
Wheelchair stabilizer for children	100.020	p. 54 i


# MOTOmed<sub>®</sub> gracile12

### Leg trainer

Ļ

Order number:

594.003



### Basic equipment (For all details, see equipment and functions page 23)

Pediatric safety foot shells (pair, incl. pedal cranks) Pediatric leg guides with calf shells (pair)

### Technical data

Dimensions (L×W×H):	$68-85 \times 55 \times 76-106$ cm
Weight:	25 kg
Maximum patient weight:	90 kg
Display size:	5,7" / 14,5 cm
Classification:	Ila according to MDD 93/42/EEC / MDR (EU) 2017/745

Mains voltage: 100-240 V~ / max. 120 VA, mains frequency: 47-63 Hz, class II equipment / Type BF

### Frequently used accessories (More accessories from page 40)

Designation	Order No.	Page
Plastic-coated pediatric safety foot shells (pair, incl. pedal cranks)	152.120	p. 44 (i)
Pediatric leg guides with plastic-coated calf shells (pair)	170.100	p. 45 (i)
Pedal radius quick adjustment for pediatric safety foot shells	588.020	p. 46 (i)
(leg trainer, pair)		
Self-operated foot holders (pair)	598.000	p. 45 (i)
Safety foot shells for adults, gracile12 (pair, incl. pedal cranks)*	586.100	p. 44 (i)
Flat stand (flat steel)	587.100	p. 55 (i)

<sup>\*</sup>Only available in combination with flat stand (order no. 587.100)



# Contact us

Email: info@neurorehabrecovery.com

Telf: 1 (866) 476 - 6897

Website: www.neurorehabrecovery.com







