


Unlock the Power
of Movement:

EXOATLET




EXOATLET EA2



Rapid don-off

Instant Footwear
Easy and fast to wear with your shoes on (Time to wear: not more than 15 sec)



Safety

Immediate stop button and handle for therapist
Excellent stability is provided by the handles during patient gait training
Also if the robot needs to be stopped, just press the stop button to halt it immediately



High-Capacity Battery

Equipped with a high-capacity battery with CE mark, the robot enables to implement robot gait training without its recharging for full daytime


- A more refined design
- More upgraded safety and convenient functions
- **FDA, Korea FDA & CE**

Indication

Stroke, TBI (Traumatic Brain Injury), SCI (Spinal Cord Injury), Parkinson Disease, Multiple Sclerosis


Simple and Quick Size Adjustment

Length of pelvis
Pull and turn the two cam levers to easily adjust the desired size



Length of thigh
Pull and release the cam lever to adjust the desired size

Length of shank
Pull and release the two cam levers to adjust the desired size




A Variety of Gait Pattern Modes

Minute bar setting of the application for robot operation in tablet PC enables patient to get customized gait training in accordance with indication and status.




Accessory
**Movable
Suspension**


This movable suspension product is an equipment used for effective assist for the purpose of expanding indications, granting therapist's degree of freedom, and enhancing stability



Green Lever on the wheel that can be fixed straight when walking forward and back
Red lever on the wheel to lock the suspension from moving



By connecting the equipment and the device, it presents the patient from falling and improves the sense of stability



Fast patient transfer is possible by using the car's seat belt method

Accessory
Crutches

Length can be adjusted according to the patient's body size. Patients can train themselves by combining smart crutch



Bluetooth



Support 4 Modes

Sitting down Standing up Walking in place Walking




Various Gait Patterns

Provides a variety of walking modes that can be set according to body size





Empowering Lives
Through Advanced
Neuro Rehabilitation



NEURO REHAB RECOVERY

Get In Touch.

Phone

1 (866) 476-6897

info@neurorehabrecovery.com

Website

www.neurorehabrecovery.com