



01

Passive Flexion & Extension:

The rehabilitation glove will drive the affected hand to grasp the small ball/water cup and other objects. Patients could intelligently adjust the time of flexion and extension.

02

Innovative Mirror Training:

With the help of syrebo™ rehabilitation glove, the healthy hand drives the affected hand to move synchronously, activates the mirror neurons, promotes the autonomic recovery of the brain and speeds up the rehabilitation process of hand function.

03

ADLs Training

The method based on motor control theory emphasizes the specificity of sports training, and reshapes the self-care ability through repeated life scene training.

Syrebo™ Soft Rehabilitation Glove is Your Best Companion!



Syrebo™ hand rehabilitation glove



Specification

Product name	Model	Size	Charge voltage	Charge Current
Soft Rehabilitation Glove	SY-HRC10	105*95*60(mm)	5V	3A

Standard accessories

Name	Qty (set)	Name	Qty (set)
Host	1	Power-assisted glove	1
Data glove	1	Adapter	1
User manual	1	Rehab ball	1
Hand pillow	1	Qualified certification	1